

## Product Data Sheet

# TruVision™ Poly-HRP IHC Detection Reagents

(anti-mouse/rabbit)

Catalog number : IHC-700-15, IHC-700-110

### Introduction

The **TruVision™** IHC Detection System is based on the site-directed clustering/coupling chemistry by which Poly-HRP conjugates are created as “*polymerized peroxidase clusters*” attached to detection antibodies. The “*polymerized peroxidase clusters*” on detection antibody offers outstanding sensitivity and signal intensity. This novel approach employed to prepare enzyme-labeled secondary antibodies avoids addition of linear polymer backbone molecules or biotin-avidin complex to reduce non-specific bindings in IHC staining. The cluster nature of polymerized peroxidase allows better permeation of the reagents to the target sites in tissue sections.

### Intended Use

For Research Use Only

### Applications

Immunohistochemistry: paraffin-embedded human tissue sections and frozen sections.

### Kit Includes:

Cat. No. IHC-700-15: The reagents provided are sufficient for 150 slides (based on 100 µl/slide)

Blocking Reagent (ready to use)	15 ml
Prediluted Secondary Reagent (ready to use)	15 ml
Anti-mouse/rabbit Poly-HRP Conjugate (ready to use)	15 ml

Cat. No. IHC-700-110: The reagents provided are sufficient for 1100 slides (based on 100 µl/slide)

Blocking Reagent (ready to use)	110 ml
Prediluted Secondary Reagent (ready to use)	110 ml
Anti-mouse/rabbit Poly-HRP Conjugate (ready to use)	110 ml

### Storage and Stability

All reagents are to be stored at 2-8 °C. The expiration date is specified on each reagent label.

### Instruction

Before beginning the staining procedure, the tissue sections should be treated via **deparaffin/rehydration** and **antigen retrieval** process. The endogenous peroxidase can be blocked with 0.3 % H<sub>2</sub>O<sub>2</sub> in TBS wash buffer. (incubate for 20 min at room temperature). Rinse 3 times with wash buffer after the blocking. Do not allow tissue sections to dry out during the staining procedure.

#### Staining Procedure:

1. Apply **Blocking Reagent** to the sections, incubate for 20 min at room temperature.
2. Drain slides to remove the blocking solution and wipe around the sections with tissue paper.
3. Apply prediluted primary antibody and incubate for 30-60 min at room temperature (or follow the instruction provided by primary antibody suppliers).
4. Rinse off the excess primary antibody with wash buffer, 3 x 2 min, with gentle agitation. Wipe around the sections with tissue paper.
5. Apply **Prediluted Secondary Reagent**, incubate for 20 min at room temperature.
6. Rinse off the reagent with wash buffer, 3 x 2 min. Wipe around the section with tissue paper.
7. Apply **Poly-HRP Reagent** to the slide, incubate for 20 min at room temperature.
8. Rinse off the excess primary antibody with wash buffer, 3 x 2 min, with gentle agitation. Wipe around the sections with tissue paper.
9. Apply DAB or AEC chromogen substrate.
10. Counterstain and mount.