



Product Data Sheet

HRP Conjugated Anti-Human IgG, Ready-to-use

Catalog number : GH9271

Product Description

This ready-to-use HRP conjugate reagent is designed for use in ELISA and Immunoassay development. It should be used with human primary antibodies in immunoassays. Prepared in HRP conjugate diluent, the concentration of HRP conjugate has been optimized for best detection performance. With stable HRP enzyme activity, this anti-human IgG HRP conjugate reagent offers convenience, consistency and reliable detection performance.

Provided As

Product No. GH9271-120: 60 ml (ready-to-use reagent)

Product No. GH9271-C: bulk qty.

Applications

ELISA

Product Use

The HRP conjugate concentration has been optimized as ready-to-use reagent for optimal detection performance.

After washing out unbound primary antibody.

Add 100 μ l into each well and incubate at room temperature for 30 minutes.

Intended Use

For Research Use and Manufacturing Purpose Only

Analysis and Test

The quality of this product is assured via ELISA, Immunoblotting and thermal stability tests.

Storage

Store at 2-8 °C.