



Product Data Sheet

Antigen Retrieval Buffer (10X), Sodium Citrate, pH 6.0

Catalog number : SU-006

Product No. SU-006-100: 100 ml

Product No. SU-006-1L: 1000 ml

Product No. SU-006-C: bulk qty.

Product Use

Dilute the concentrate to 1X as working solution. Proceed with normal target unmasking step using this 1X antigen retrieval buffer.

Product Description

The antigen retrieval buffer contains sodium citrate, pH 6, and is supplied as a concentrate format (10x). It is commonly used for unmasking targets of FFPE tissue sections.

Intended Use

For Research Use Only

Storage

Store at 2-8 °C.