



## Product Data Sheet

### Antigen Retrieval Buffer (10X), EDTA, pH 8.0

Catalog number : SU-007

---

Product No. SU-007-100: 100 ml  
Product No. SU-007-1L: 1000 ml  
Product No. SU-007-C: bulk qty.

#### Product Use

Dilute the concentrate to 1X as working solution. Proceed with normal target unmasking step using this 1X antigen retrieval buffer.

#### Product Description

The antigen retrieval buffer, EDTA solution, pH 8 is commonly used in unmasking targets in IHC procedure.

#### Intended Use

For Research Use Only

#### Storage

Store at 2-8 °C.