



Product Data Sheet

Antigen Retrieval Buffer (10X), Tris-EDTA, pH 9.0

Catalog number : SU-008

Product No. SU-008-100: 100 ml
Product No. SU-008-1L: 1000 ml
Product No. SU-008-C: bulk qty.

Product Use

Dilute the concentrate to 1X as working solution. Proceed with normal target unmasking step using this 1X antigen retrieval buffer.

Product Description

Supplied as concentrate format (10x), Tris-EDTA, pH 9, is commonly used in unmasking targets in IHC procedures.

Intended Use

For Research Use Only

Storage

Store at 2-8 °C.