

Product Data Sheet

Antibody Diluent

Catalog number : SU-003-15, SU-003-110

Product Description

The ready-to-use antibody diluent is recommended for use in Immunohistochemistry, ELISA immunodetection assays, for diluting primary antibodies, HRP conjugates or other secondary antibodies. This product may also be used as initial blocking reagent for initial block steps in IHC or ELISA procedures.

Cat. No. SU-003-15: 15 ml ready-to-use diluent solution.

Cat. No. SU-003-110: 110 ml ready-to-use diluent solution.

Intended Use

For Research Use and Manufacturing Purpose Only

Storage and Stability

Store at 2-8 °C, stable for at least 12 months from the date received.

Product Use

Antibody diluent is recommended for use in diluting primary and secondary antibody in IHC, ELISA or Western Blotting procedures.

Gently mix concentrated antibody in diluent, avoid foam forming.