

Product Data Sheet

Anti-Mouse Poly-HRP IHC Reagent

Catalog number : IHC-5281-15, IHC-5281-50, IHC-5281-C

Product Description

This ready-to-use Poly-HRP conjugate reagent is designed for use in immunohistochemistry procedure. It should be used with mouse primary antibodies. Prepared in HRP conjugate diluent, the concentration of Poly-HRP conjugate has been optimized for best detection performance. With stable HRP enzyme activity, this anti-mouse IgG Poly-HRP conjugate reagent offers convenience, consistency and reliable detection performance.

<u>Cat. No. IHC-5281-15:</u>	15 ml
<u>Cat. No. IHC-5281-60:</u>	60 ml
<u>Cat. No. IHC-5281-C:</u>	bulk qty.

Intended Use

For Research Use or Manufacturing Purpose Only

Analysis and Test

The quality of this product is assured via IHC test

Storage and Stability

Store at 2-8 °C, stable for at least 12 months from the date received.

Applications

Immunohistochemistry
ELISA
Western Blotting

Product Use

The HRP conjugate concentration has been optimized as ready-to-use reagent for optimal detection performance.