

Product Data Sheet

Goat Anti-Mouse IgG (H+L), Affinity-Purified

Catalog number : GM5200-1, GM5200-10, GM5200-C

Product Description

This affinity-purified antibody binds to heavy and light chains on mouse IgG. It reacts with all mouse IgG subclasses. Supplied as a concentrated preparation in PBS buffer, this goat anti-mouse IgG (H+L) antibody does not contain preservatives or any carrier protein, and is available in large volumes for further manufacturing use.

Intended Use

For Research Use or Manufacturing Purpose Only

Storage

Store at -20 °C

Provided As

<u>Cat. No. GM5200-1:</u>	1 mg in PBS.
<u>Cat. No. GM5200-10:</u>	10 mg in PBS.
<u>Cat. No. GM5200-C:</u>	bulk qty.

Applications

Western Blotting
ELISA
Immunohistochemistry

Analysis and Test

The quality of this product is assured via ELISA and Immunoblotting tests.