

Product Data Sheet

Goat Anti-Rabbit IgG (H+L), HRP Conjugated

Catalog number : GR2201-5, GR2201-C

Product Description

This HRP-conjugated goat anti-rabbit IgG (H+L) secondary antibody reacts with heavy and light chains on rabbit IgG. The HRP conjugate is prepared via our proprietary production process, and offers high sensitivity and clean detection performance. Supplied as liquid in screw cap-vials, this HRP conjugate is stable at 4 °C.

This secondary antibody is recommended for use in routine immunodetection assays. Working dilutions should be determined by the investigators to obtain the optimal detection performance. For example, working dilution may begin with 1:3000 for Western Blotting, 1:5000 for ELISA.

Cat. No. GR2201-5: 0.5 ml
Cat. No. GR2201-C: bulk volume

Intended Use

For Research Use or Manufacturing Purpose Only

Storage and Stability

Store at 2-8 °C, stable for at least 12 months from the date received.

Analysis and Test

The quality of this product is assured via ELISA and Immunoblotting tests.

Applications

Western Blotting
ELISA

Recommended Working Dilutions

Optimal working dilutions should be determined by the investigators to obtain the optimal detection performance.

-Western Blotting: 1:3000
-ELISA: 1:5000